



Safe Surgery any time, any place®

Presenting
the award-winning
SurgiField™ System



SurgiBox Inc. is a Delaware registered C-Corp

NATO and the U.S. military intend to push lifesaving surgical capability far forward into austere and contested environments.

Future conflicts will **stretch evacuation** timelines. NATO and the U.S. Military recognize that the golden hour cannot be guaranteed. What can be done ?

- Push surgical capability forward: Both NATO and U.S. doctrine emphasize moving damage-control surgery to Role-2/FST-level echelons in austere environments when evacuation is delayed or contested (LSCO / A2AD contexts).
- Prolonged Field Care (PFC): JTS Clinical Practice Guidelines mandate forward interventions to sustain casualties beyond the "golden hour" when MEDEVAC is unavailable.
- Mobile surgical/resuscitative teams: Recent NATO/US papers highlight the need for small, transportable Role-2 augmentations capable of parallel surgeries in distributed operations.

The SurgiField™: A **rapidly** deployed, **cost-effective** and **ultra-portable** surgical environment designed with US **Army** and **Air Force** surgeons:

Torso




Limb



- protects the patient's surgical site from the outside environment and medical personnel;
 - the medical provider from the patient's bodily fluids;
-
- Role-2 augmentation in a backpack
 - Rapid setup for parallel surgeries
 - Enables PFC & en-route resuscitation
 - Low logistics footprint, NATO interoperable
 - Force multiplier for SOF, MEU, Littoral ops



The SurgiField™



Weight
< 6 kgs



Inflates
~ 120s



Each charge
4+ hours



Dual Patient-Provider
Protection

"When you're operating with artillery fire in the distance and debris falling around you, you learn that medicine isn't about perfect conditions - it's about creating a moment of safety in chaos. **The SurgiField gave us something I never thought possible in that building:** better sterile conditions. For those 30 minutes, with my hands working inside that protective barrier, I could focus entirely on saving this young man's life instead of worrying about contamination. The system was intuitive even under pressure - clear visibility, easy access, and most importantly, it created a clean surgical field where none should have existed.

In conflict medicine, we're always improvising, but this felt like having a real operating room materialize around my patient. That's the difference between losing someone to infection days later and giving them a fighting chance to heal."

Burmese Army Surgeon,



Shrapnel removal in Burma

Real Life Impact In Ukraine through
Emergency Use Authorization: "This
equipment helps to save, firstly, time for
processing the operating field, and
secondly, it helps to ensure sterility in
aggressive environments. Namely, in
basements, tents, under the conditions of
transportation by sanitary transport. This is
what provides **immediate operative access**
to a serious injury, ensures the preservation
of the life of such an injured person where
it is not possible to provide qualified
assistance in the operating room."
(translated from original Ukrainian),

Lieutenant General, Ukraine



[Check out the full video here !](#)



До військового шпиталю, що виконує завдання у Північній операційній зоні, передали портативні інноваційні системи «SurgiBox» для невідкладних хірургічних операцій в польових умовах. Про це [повідомив](#) командувач Об'єднаних Сил ЗС України генерал-лейтенант Сергій Гасій.

"If you're going to go to a surgically underserved country, you're going to try to find a hospital somewhere and that hospital may or may not have the infrastructure to support clean surgery. But when you think about the SurgiField, suddenly you don't need any of that. You can go to the patient instead of the patient coming to you."



*Prof. David King, Trauma Surgeon,
Associate Professor, Massachusetts General
Hospital and US Army Surgeon*



Japanese Surgeon and Military Simulating the SurgiField

War Settings



The SurgiField was designed with trauma surgeons from US Army and US Air Force



Field Surgery

- Field Hospitals / Naval Ships: Parallel surgeries, higher throughput.
- Prolonged Field Care: On-site surgery when evacuation impossible.
- Littoral Combat Ships: Makeshift operating rooms, major upgrade for surgeons.



Evacuation

- CASEVAC: Manage complications in flight, protect patients/crew.
- Medevac: Stabilization + transport, wound protection

Emergency, Austere and Remote Setting Scenarios



Emergency Settings (e.g. conflict and disaster zones) to enable **the transformation of any space** into a safe surgical environment with multiple “pop-up” operating theaters including:

- Existing infrastructure such as medical facilities - even those that have been damaged
- Infrastructure that can be set up on the go such as in tents, abandoned buildings and sheds



Austere and Remote Settings (e.g. refugee camps, rural areas, mines, cruise ships) to enable the deployment of surgical missions in:

- Terrain (e.g. mountains and mines) where transporting patients is difficult
- Remote areas where no or poor operating theater infrastructure exists and/or reliance on distant referral hospitals is impractical

The SurgiField™ has been developed with scientific rigor with results published in peer-reviewed journals

MILITARY MEDICINE, Vol. 186, January/February Supplement 2021

A Paradigm Shift in Critical Care Infrastructure in Complex Settings: Evaluating an Ultraportable Operating Room to Improve Field Surgical Safety

Debbie L. Teodorescu, MD^{1,1,1,1}; Stephen Chakima, MS^{1,1}; CPT Asad Moten, MC, USA^{1,1,1,1,1,1}; Mike H.M. Teodorescu, DDA^{1,1,1}; Major [F] Hoch, MD^{1,1,1}; Samantha Gallimore-Briggs, BS^{1,1,1}; Karim Maitz, MD, MS^{1,1}; Capt Robert J. Crowley, MC, USA^{1,1}; COL David R. King, MC, USA^{1,1,1,1,1}

ABSTRACT

Introduction:

Scarcity of operating rooms and personal protective equipment in forward field settings make surgical infections a potential concern for combat mortality and morbidity. Surgical and transport personnel also face infectious risks from bodily fluid exposures. Our study aimed to describe the rapid, proof-of-concept testing of the SurgiField technology as an inflatable, mobile environment that addresses the aforementioned problems, fits on generators and backpacks, and disposes over incinerators.

Materials and Methods:

The SurgiField environmental control unit and inflatable enclosure were optimized over five generations based on negative feedback from stakeholders experienced in surgery in austere settings. The airflow system was developed by analytic modeling, verified through in vitro modeling in a 54 (5000) and confirmed with primary smokestack checking. Particulate counts evaluated the enclosure's ability to contain and mitigate most exposures to potentially infectious contaminants from the surgical field in various settings. SurgiField enclosures were setup over a mannequin's torso in a configuration and position for either thoracic or abdominal surgery. A particle counter was centrally positioned in maneuver and deployment positions, as well as lateral field positions. This setup was repeated with open ports exposing the enclosure to the external environment. To simulate stress scenarios, sampling was repeated with enclosure measurements during an increase in external particulate concentration.

Results:

The airflow technology effectively kept contaminants away from the incision and maintained a pressure differential to reduce particulate entry. Flowing airways demonstrated that even when ports were opened in the external environment, high containment levels, the enclosed surgical field consistently captured all particles, even at all positions. Time from key opening to incision averaged 54.5 seconds, with the rate-limiting step being connecting the environmental control unit to the incision. The portable kit weighed 5.9 kg.

Conclusions:

Analysis, in silico, and mechanical airflow modeling and benchtop testing have helped to quantify the SurgiField system's capability to contain and maintain an operating room quality surgical field within the enclosure as well as protecting the surgical team outside the enclosure. More recent and ongoing work has focused on specifying optimal air settings in the monthly check of air, expanding support for streamlined procedures, automating airflow control, and accelerating system setup. SurgiField's ultimate goal is to take timely, safe surgery to patients in even the most austere of settings.

¹Department of Research and Development, SurgiBox Inc., Cambridge, MA 02142, USA.
²Department of Medicine, Massachusetts General Hospital, Harvard Medical School, Boston, MA 02114, USA.
³USLA, Massachusetts Institute of Technology, Cambridge, MA 02139, USA.
⁴Department of Defense, National Capital Consortium, Defense Health Agency, Washington DC 20305, USA.
⁵Department of Lab, Health Research International, Houston, TX 77060, USA.
⁶Department of Information Systems, Boston College, Chestnut Hill, MA 02465, USA.
⁷United States Military Medical Center, Loma, CA 92556, Northridge.
⁸Medical Service, United States Air Force, Washington, DC 20330, USA.

MILITARY MEDICINE, Vol. 186, January/February Supplement 2021

295

Early-stage innovation report

OPEN ACCESS

User-centred innovation to advance safe surgery at the point of need: report on development of the SurgiField system

Debbie Lin Teodorescu¹, U¹ Soththar Jongsakulapadda¹, Atif Rakin¹, Mike H M Teodorescu¹, Macaulay Kenney¹, Julian Greene¹, Steven Bokhan¹, Hajar Mahdian¹, Asad Moten¹, Omid Abo¹, Robert D Smalley¹, Ander Dorken Callagott¹, Arnav Mahajan¹, Robert D Smalley¹, Daniel Frey¹, David R King¹, Meena N Chertan¹

INTRODUCTION

In conditions ranging from pneumothorax to acute compartment syndrome to fractures to traumatic and degenerative joint disease to appendicitis to obstructed labour to postpartum haemorrhage, timely surgical intervention is the deciding factor among emergency mortality or long-term disability. Consequently, access to quality surgical care is a global priority across the resource spectrum. To date, this, an estimated 18 million people die annually from conditions that safe surgery could address. Even more experience permanent disability from infirmed, inaccessible, or unsafe surgery: generating a gross domestic product loss of US\$20 billion from lost annual production¹. For those requiring reconstructive, sophisticated surgeries, the vast majority of surgical needs globally are captured by the simple, scalable, full-spectrum surgical procedures.

Surgical care has traditionally relied on consistent availability of three resource categories: (1) functional equipment and materials, (2) appropriate facility space and (3) medical staff^{2,3}. Given the practical challenges of disseminating these resources in many settings, including but not limited to conflict settings, disaster zones, low-to-middle-income countries and rural regions, surgical care availability is often limited to centralized referral settings⁴. Even higher-resource settings regard cost, complications and quality improvement as ongoing challenges to care equity and sustainability. For reference, in some operating theatres cost an average of US\$27 per minute to run⁵. High equipment and facility costs are driven in part by maintaining the sterility

WHAT ARE THE NEW FINDINGS

→ Proven, iterative development of individual elements of a portable technology may make legal as well as highly optimized innovation with innovation, ongoing care feedback, however, emerging innovation risk, such as the complex integration, is critical to overall success of the innovation.

→ The SurgiField system is able to provide both passive and active barriers against environmental contamination during a broad range of surgical procedures, particularly for so-called 'bulletproof' procedures that are used intensively in forward theatre areas to surgical care.

→ A flexible, sealed, full-sized enclosure could be made efficient to deploy and take down, transport and equipped with typical surgical field access requirements. This could be implemented with a conveniently sized package similar to large surgical drapes and readily for any major surgical procedure.

Check for updates

Submitted to BMJ Innovations 2020; accepted 2020; published 2021

© 2021 Teodorescu et al. Published by Oxford University Press on behalf of the Association of Military Surgeons of the United States 2021. This work is written by US Air Force Government employees and is in the public domain in the US.

Teodorescu DL, et al. BMJ Innovations 2021;17(1):1-11. doi:10.1093/bjino/ibz008

317

European Journal of Trauma and Emergency Surgery

https://doi.org/10.1007/s00120-020-02470-w

ORIGINAL ARTICLE

Moving toward point-of-care surgery in Ukraine: testing an ultra-portable operating room in an active war zone

Emanuele Lagazzi^{1,2,3}, Debbie Lin Teodorescu^{1,4}, Elias Argandhyer⁵, Michael Alexander Samotkovski¹, David Richard King^{1,6}

Received: 26 September 2021 / Accepted: 19 November 2021
© The Author(s), under exclusive license to Springer Verlag GmbH Germany 2021

Abstract

Purpose: In conflict zones, providers may have to decide between delaying time-sensitive surgeries or performing operative interventions in the field, potentially subjecting patients to significant infection risks. We conducted a single-arm crossover study to assess the feasibility of using an ultraportable operating room (U-OR) for surgical procedures on a porcine cadaver abdominal traumatic injury model in an active war zone.

Methods: We enrolled participants from an ASSET-type course designed to train Ukrainian surgeons before deployment to active conflict zones. They performed three standardized consecutive abdominal surgical procedures (liver, kidney, and small bowel injury repair) with and without the U-OR. Primary outcomes included surgical procedure completion rate, procedure time, and airborne particle count at the start of surgery. Secondary survey-based outcomes assessed surgery task load index (SURT-X) and perceived operative factors.

Results: Fourteen surgeons performed 76 surgical procedures (38 with the U-OR, 38 without the U-OR). The completion rate for each surgical procedure was 100% in both groups. While the procedure time for the liver injury repair did not differ significantly between the two groups, the use of the U-OR was associated with a longer time for kidney (155 vs. 96 s, $p=0.002$), and small bowel (229 vs. 101 s, $p=0.004$) injury repair. The average airborne particle count within the U-OR was substantially lower compared to outside the U-OR (6,753.85 vs. 237,282 and, $p<0.001$). There was no statistically significant difference in SURT-X for procedures performed with and without the U-OR.

Conclusion: The use of the U-OR did not affect the procedure completion rate or SURT-X. However, there was a marked difference in airborne particle counts between inside and outside the U-OR during surgery. These preliminary findings indicate the potential feasibility of using a U-OR to perform abdominal damage-control surgical procedures in austere settings.

Keywords: Trauma · Emergency general surgery · Austere environment · Disaster medicine · Ukraine

Michael A. Samotkovski and David R. King are co-senior authors with equal contribution.

✉ David Richard King

dking@bgh.harvard.edu

Emanuele Lagazzi

elagazzi@bgh.harvard.edu

Debbie Lin Teodorescu

dteodorescu@bgh.harvard.edu

Elias Argandhyer

eliasargandhyer@bgh.harvard.edu

Michael Alexander Samotkovski

msamotkovski@bgh.harvard.edu

Published online: 04 January 2021

Springer

- [Military Medicine \(Link Here\)](#)
- Results of the benchtop testing

- [BMJ Innovations \(Link Here\)](#)
- Results of the co-development process with surgeons and medical personnel

- [Pubmed \(Link Here\)](#)
- Results of a study in Lviv Ukraine

Shipping Information

	SKU	Package Weight (cms)	Package Dimensions	Shelf Life/Warranty
	SB-T1-02-CE Carton of CE SurgiBubbles (6x)	62.35 x 51.4x35.3	8kgs	2 Years
	EB-T1-02-EU: 1 Smart control Module 1 Smart Charger EU 2 Lithium Batteries 2 Pop-Up Frames	Box 1: 45.0 x 32.2 x 28.6 Box 2: 97 x 10 x 10	Box 1: 6.75 kgs Box 2: 0.3 kgs	2 Years



TIME Magazine's Best Inventions
of 2023



Safe Surgery Any Time, Any Place®

sales@surgibox.com